


26-27
APRIL
2018

2nd EuroTauMeeting

Or how to know everything about Tau protein and associated diseases ?

LILLIAD Learning Center Innovation
Avenue Henri Poincaré - Villeneuve d'Ascq

Mirita



ABSTRACTS & PROGRAM BOOK

Organized by

Scientific committee

Luc BUEE, Amrit MUDHER, Jean-Pierre BRION

Miguel MEDINA and Jesus AVILA

Local committee

Sophie HALLIEZ

Valerie BUEE-SCHERRER

Charlotte CAILLOT

and the Alzheimer & Tauopathies Team









Lille



Conference secretariat
charlotte.caillot@inserm.fr



Metro line 1 (yellow)
Direction 4 cantons - Stade P. Mauroy
Station «cité scientifique»



P08: Inhibition of protein phosphatases by okadaic acid induces expression of high molecular weight phospho-tau-immunoreactive protein species in neuroblastoma SH-SY5Y cells.

Mirta Boban, Terezija Miškić, Mirjana Babić Leko and Goran Šimić.

Department of Neuroscience, Croatian Institute for Brain Research, University of Zagreb Medical School, Zagreb, Croatia.

e-mail: mirta.boban@hiim.hr

A key feature of Alzheimer's disease (AD) is aggregation of microtubule-associated protein tau in the neurofibrillary tangles (NFT) in the brain. NFT tau is characterized by abnormally high phosphorylation, which may result from the upregulated activity of protein kinases and downregulation of protein phosphatases. To investigate tau under the condition of protein phosphatase impairment, we treated neuroblastoma SH-SY5Y cells with okadaic acid (OA), an inhibitor of protein phosphatases and analyzed total cell lysates with phospho-tau and total tau antibodies using immunoblot. In addition to the well-described 50-65 kDa tau isoforms, we observed that both undifferentiated and retinoic acid- and brain derived neurotropic factor-differentiated SH-SY5Y cells treated with OA express high molecular weight protein species immunoreactive with anti-tau-pS202 and -pS396 antibodies. The apparent molecular weight of 100 kDa indicated a possibility of tau dimer. In support, high molecular weight tau immunoreactive proteins (HMW-TIP) were detected in a heat-stable fraction. However, we were unable to detect HMW-TIP using anti-total tau antibodies. This could be due to protein truncation or epitope masking within the oligomer, or a possibility that HMW-TIP represents a tau-unrelated protein. Our biochemical characterization showed that HMW-TIP were stable under reducing conditions and in the presence of strong denaturing agents, such as urea and guanidine, as well as upon alkaline phosphatase treatment. In conclusion, we show that protein phosphatase inhibition by okadaic acid induces the appearance of HMW-TIP, which may represent tau oligomer or tau cross-reactive phospho-proteins.